



# BANTAM | PRO™ -G

## *A952-G Electrosurgical Generator Package*



### **COMPLETE**

The BOVIE BANTAM | PRO A952-G is packaged with the essential accessories needed to perform an OB/GYN LEEP procedure.

### **CONVENIENCE**

The full system is situated on a mobile stand for effortless portability.

### **AFFORDABLE**

The BOVIE A952-G is an extremely affordable electrosurgical generator package with smoke evacuation.

### **LLETZ LOOP ELECTRODES**

A variety of electrodes, plus LLETZ Loop Electrodes, are included.

### **RELIABILITY**

The BOVIE A952 unit is backed by Bovie's four (4) year manufacturer's limited warranty.

- High Frequency Desiccator with 50 watts of Cutting capabilities
- Two machines in one convenient package
- Adjustable in 1/10th watt increments in the first 10 watts, for the most accurate output available
- No grounding pad required when working in the desiccation mode
- When in electrosurgical generator mode, the unit features Bovie NEM™ (Neutral Electrode Monitoring) pad sensing technology which monitors the return electrode for optimum patient protection.



To learn more visit [BovieMedical.com/bantampro-g](http://BovieMedical.com/bantampro-g)

# BOVIE® BANTAM | PRO™ - G ELECTROSURGICAL GENERATOR PACKAGE


## EQUIPMENT INCLUDED



A952



A902



ES01, ES02, ES20



A804, A805, A806, A807, A806DE, A807DE



A1252C



ESRE



A910ST



A910



A812-C



ESLK



A837



SE02



ES07



786T



SERF

### Style # A952-G

#### Standard Equipment

- A952** (1) electrocautery generator,
- A902** (1) handpiece
- A804** (5) non-sterile sharp dermal tips
- A805** (2) sterile sharp dermal tips
- A806** (5) non-sterile blunt dermal tips
- A807** (2) sterile blunt dermal tips
- A806DE** (2) non-sterile DERM-Elite™ Premium blunt dermal tips
- A807DE** (2) sterile DERM-Elite Premium blunt dermal tips
- ES01** (3) blade sterile disposable electrodes
- ES02** (1) needle sterile disposable electrodes

- ES20** (1) ball sterile disposable electrodes
- A910ST** (2) sterile disposable handpiece sheaths
- A910** (2) non-sterile disposable handpiece sheaths
- A1252C** (2) reusable grounding cord
- ESRE** (5) disposable split grounding pads
- A837** (1) wall mount kit
- 110 VAC Hospital-grade power cord\*
- User's guide (MC-55-238-001)\*

#### • Smoke Shark™ II system

- + **SE02** Smoke Shark with pneumatic footswitch
- + **SF35** 35-hour filter
- + **786T** 7/8" (22.23mm) tube 6' (182.9cm) long
- + **SERF** reducer fitting\* (attaches to vaginal speculum)

#### • A812-C electrocautery mobile stand

- + **A812** mobile stand
- + **A812-E** height extender
- + **A808-T** instrument tray
- + **A812-BT** bottom tray

#### • ESLK electrocautery loop kit

- + **ES09** 10mm x 10mm loop
- + **ES11** 15mm x 10mm loop
- + **ES12** 20mm x 8mm loop
- + **ES13** 20mm x 15mm loop
- + **ES16** 10mm x 8mm square
- **ES07** five extended ball electrodes

\* Not Pictured