



# SIGNATURE

## Exam Table Series

Now Better Than Ever



EMPOWERING  
**HEALTHY**  
OUTCOMES

Improving Quality, Safety & Efficiency

# Signature Series

## Exam Tables

# Compare Models



**PROfrontstep**

5240

**PROsidestep**

5240-145

**ULTIMA**

5250

### KEY FEATURES

500lb capacity (226.8kg)	✓	✓	✓
Removable seamless upholstered cushion	✓	✓	✓
Easy-clean powder coat enamel surfaces	✓	✓	✓
Smooth moving articulating back adjustment	✓	✓	✓
Stirrups	✓	✓	✓
Sturdy steel extension leg pullout with seamless, ultra-comfort leg cushion	✓	✓	✓
Removable 12.8"x 10.3" plastic drip pan (32.5 cm x 26.1 cm)	✓	✓	✓
Pull-out all steel non-slip front step 7" high (17.8 cm)	✓	—	✓
Pull out steel non-slip side step 7" high (17.8 cm)	✓	✓	✓
21" Paper roll holder (53.3 cm)	✓	✓	✓
2 front drawers on full extension roller bearing glides 21.00" x 17.63" x 5" (53.34 cm x 44.70 cm x 12.7 cm) each	✓	✓	✓
2 reversible side storage drawers on full extension roller bearing glides 21.00" x 17.63" x 5" (53.34 cm x 44.70 cm x 12.7 cm) each	✓	✓	—
3 reversible side storage drawers on full extension roller bearing glides 21.00" x 17.63" x 5" (53.34 cm x 44.70 cm x 12.7 cm) each	—	—	✓
5 year limited warranty	✓	✓	✓

### ACCESSORIES

010 Leg extension booster pad	810 Patient Safety Straps (Factory installed only)
037 Armboard assembly (Specify patient right or left side / factory installed)	820 Stainless steel drain pan (Factory installed)
045 Accessory rail (Pair / specify head or seat section / factory installed)	822 Stainless steel drain pan with spout (Factory installed)
251 Bierhoff knee crutches	849 Drawer warmer 220V (Factory installed only)
350 Pelvic tilt	851 Drawer warmer 110V (Factory installed only)
365 I.V. pole (Factory installed)	852 Hospital grade receptacle 110V with child safe cover (Factory installed / patient left side)
375 Safety rails assembly (Pair / factory installed)	855 Hospital grade receptacle 220V with child safe cover (Factory installed / patient left side)
463 Mounting Plate for Welch Allen Light # 48955	5250F Fire retardant foam
701 Removable Seamless Upholstered Top with Antimicrobial Guard	5250B Base only
701F Removable Seamless Upholstered Cushion with CAL-133 foam	SilverSan antimicrobial powder coated base

### TABLE SPECIFICATIONS

<b>Height:</b> 31.8" (80.8cm)
<b>Length:</b> 58.5" - 77.5" (148.6cm - 196.9cm)
<b>Width:</b> 28" (71.12cm)
<b>Stirrups:</b> 4 lateral positions & adjustable length
<b>Construction:</b> Galvanized steel base, 16 gauge, powder coated

### COLORS



# Features

**500 lb. patient weight capacity**

**Infection Control** – seamless upholstery, removable for cleaning and replacing

**Seamless Surfaces** – Seamless drawers and easy-clean enamel surfaces

**Durability & Reliability** – Solid Steel Construction

**Efficiency** – Up to 5 Large Storage Drawers

**7" step height** for stability

# Praise

“After thousands of patient visits in our busy primary care practice, the UMF Medical exam table withstands the test of time. You don’t find this quality anywhere else.”

---

“These tables are built to last”

---

“This new table looks so inviting, sturdy and comfortable! And a big guy like me doesn’t feel timid to get on it.”



# The Signature Difference



Next Generation Comfort Cushion for patient comfort

Seamless, removable cushion top for easy cleaning

Largest drawer capacity in industry


Reversible drawers

Stirrups are infinitely adjustable for maximum comfort

500 lb (226 kg) patient weight capacity

Slim-Line paper roll holder

New smooth moving  
pneumatic backrest



**Four or five drawers**

**Side or front step**



**Infinite adjustability**



The **5200 Signature Series Exam Table** is designed to optimize the patient and caregiver experience. A perfect blend of style and function, this patient-centric table is a mainstay in busy primary care practices. Your patient population is now changing, and every detail matters. **That's why UMF Medical is setting a new benchmark for infection control, smart efficiencies, and lasting value.** By asking what's going to improve your patients' experience, UMF Medical has arrived at an exam room solution that keeps you taking care of patients, day after day.

## Empowering Healthy Outcomes

### **UMFmedical**

1316 Eisenhower Boulevard  
Johnstown, PA 15904-3393  
Phone: 814-266-8726  
Toll Free: 800-638-5322  
Fax: 814-266-1870  
[www.UMFmedical.com](http://www.UMFmedical.com)

